

Eppendorf Advantage

Advance your research through insights,
solutions, optimal outcomes, and
special offers.

JANUARY 1 – JUNE 30, 2025





Eppendorf Advantage Promotions

Advantage means putting yourself in a better position, and that's exactly what we want for you.

The work you do is more important now than ever before, and the samples you produce more precious. You have been able to rely on Eppendorf to help deliver accurate, reliable results, but we want to be more than your provider of consumables, instruments, and services.

We are also here to be your partner and resource center. Our bi-annual publication includes articles and promotional offers. We hope to help give you the Advantage to succeed in the lab!

Ordering information

Description

Part Num.

Thermocyclers Save 25% on Mastercyclers®

1 x Mastercycler® X50s, 110-230V 96-well thermoblock, silver with control panel EU and US plug and 2 x Mastercycler® X50i, 110-230V 96-well thermoblock, silver without control panel EU and US plug	EP87-911
------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------

1 x Mastercycler® X50a, 110-230V 96-well thermoblock, aluminum with control panel EU and US plug and 2 x Mastercycler® X50l, 110-230V 96-well thermoblock, aluminum without control panel EU and US plug	EP87-912
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------

1 x Mastercycler® X50h, 110-230V 384-well thermoblock, aluminum high pressure lid with control panel EU and US plug and 2 x Mastercycler® X50t, 110-230V 384-well thermoblock, aluminum high pressure lid without control panel EU and US plug	EP87-913
------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------

Plates & Accessories Save 50% on Eppendorf twin.tec® *real-time* PCR plates, and 15% on Masterclear® *real-time* PCR film

2 x twin.tec® <i>real-time</i> PCR Plate 96, unskirted, low profile, white, 20 pcs.	EP87-907
-------------------------------------------------------------------------------------	----------

2 x twin.tec® <i>real-time</i> PCR Plate 96, semi-skirted, high profile, blue, 25 pcs.	EP87-908
----------------------------------------------------------------------------------------	----------

2 x twin.tec® <i>real-time</i> PCR Plate 96, skirted, low profile, white, 25 pcs.	EP87-909
-----------------------------------------------------------------------------------	----------

Masterclear® <i>real-time</i> PCR Film (self-adhesive), 100 pieces	EP87-910
--------------------------------------------------------------------	----------

Centrifuge Bundle Save over 30%

Centrifuge 5920 R, keypad, refrigerated, 120 V/50 – 60 Hz (US), includes Rotor S-4x Universal-Large, 2 adapters, 2 caps	EP87-916
-------------------------------------------------------------------------------------------------------------------------	----------

Centrifuge Concentrator Bundle Save over 30% on the Centrifuge Concentrator Bundle

Vacufuge® plus system with integrated diaphragm vacuum pump, incl. integrated pump, 120 V/50 – 60 Hz (US), w/ Rotor F-45-48-11 & Rotor F-35-6-30; Eppendorf Conical Tubes SnapTec® 50	EP87-914
---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------

Thermomixer Save over 16% on the Eppendorf Thermomixer®

Eppendorf ThermoMixer® C, 100 – 130 V/50 – 60 Hz; Eppendorf SmartBlock 1.5 mL, incl. Transfer Rack 1.5/2.0 mL	EP87-918
---------------------------------------------------------------------------------------------------------------	----------

Ordering information

Description	Part Num.
Centrifuge & Pipette Bundles Save over 24% on the Microcentrifuge 5425 Bundle, and over 39% on a Centrifuge 5430 Bundle	
Microcentrifuge 5425, non-refrigerated, keypad control panel, 120 V/50 – 60 Hz, Rotor FA-24x2; includes 3x Eppendorf Research® plus (one-channel) pipettes with 3x epT.I.P.S.® Box 2.0 (x96 pipette tips)	EP87-896
Microcentrifuge 5425, non-refrigerated, knob control, 120V/50 – 60 Hz, Rotor FA-24x2; includes 3x Eppendorf Research® plus (one-channel) pipettes with 3x epTIPS Box 2.0 (x96 pipette tips)	EP87-900
Centrifuge 5430 R, rotary knobs, refrigerated, 120 V/50 – 60 Hz (US); aerosol-tight, incl. rotor lid; incl. wind shield and rotor lid, with 2 buckets; Eppendorf Xplorer®, 8-channel, 15 – 300 µL, orange, 100 – 240 V ±10 %, 50 – 60 Hz	EP87-915
Centrifuge 5430 R, keypad, refrigerated, 120 V/50 – 60 Hz with Rotor FA-45-30-11 aerosol-tight rotor, incl. rotor lid; Rotor A-2-MTP incl. wind shield and rotor lid, w/ 2 buckets; Eppendorf Xplorer®, 8-channel, 15 – 300 µL, orange, 100 – 240 V ±10 %, 50 – 60 Hz	EP87-919
Liquid Handling Bundles Save over 19% on Liquid Handling Bundles	
LH Essentials Electronic Package: Buy 3 Eppendorf Xplorer® single-channel pipettes (0.5 – 10 µL, 5 – 100 µL, 50 – 1,000 µL) and receive a charging carousel, discounts on a Repeater® E3 electronic dispenser, an Easypet® 3 pipette controller, and a free Pipette Manager software for efficient data management	EP87-888
LH Essentials Package 2: Features 4 Research® plus single-channel pipettes (0.1 – 2.5 µL, 2 – 20 µL, 20 – 200 µL, 100 – 1,000 µL), 3x epT.I.P.S.® Box 2.0 (each containing 96 pipette tips), a ballpoint pen, a carousel stand, an Easypet® 3 pipette controller, and a Repeater® E3x electronic dispenser with holder.	EP87-887
Manual & Electronic Pipettes Save over 25% on Manual & Electronic Pipettes	
Xplorer® plus, 8-channel, 5 – 100 µL, yellow, 100 – 240 V ±10 %, 50 – 60 Hz; Xplorer®, 8-channel, 15 – 300 µL, orange, 100 – 240 V ±10 %, 50 – 60 Hz; Xplorer® plus, 8-channel, 50 – 1,200 µL, green, 100 – 240 V ±10 %, 50 – 60 Hz	EP87-901
Eppendorf Research® plus, 1-channel, 20 – 200 µL, yellow, variable volume; Eppendorf Research® plus, 1-channel, 0.1 – 2.5 µL, dark gray, variable volume; Eppendorf Research® plus, 1-channel, 100 – 1,000 µL, blue, variable volume; Eppendorf Research® plus, 1-channel, 2 – 20 µL, yellow, variable volume incl. epT.I.P.S.® Box 2.0 with 96 pipette tips; Repeater® M4 Starter Kit, 1-channel, 1 µL – 10 mL, incl. Combitips® advanced Box 2.0, Combitips® advanced assortment pack (1 tip of each size) and Pipette Holder 2	EP87-902
Freezers Save 15% on the Innova® 42R	
Innova 42R with 3/4" orbit, 120V/60Hz; Stacking kit; Universal Platform 18" x 18"	EP87-903
Innova 42R with 1" orbit, 120V/60Hz; Stacking kit; Universal Platform 18" x 18"	EP87-904
Incubator Shakers Save on the Innova® S44i Incubator Shaker	
Innova S44i refrigerated, stackable up to 3-fold, 120 V ± 10 %/50/60 Hz (US), 1" orbit, w/ Universal platform	EP87-885
Innova S44i refrigerated, stackable up to 3-fold, 120 V ±10 %/50/60 Hz (US), 2" orbit, w/ Universal platform	EP87-886
3X Innova® S44i refrigerated, stackable up to 3-fold, 120 V ± 10 %/50/60 Hz (US), orbit 2.5 cm (1 in); 3x Universal platform; Triple stacking kit	EP87-917
CO₂ Incubators Stack & Save with the CellXpert® C170i	
2 x CellXpert® C170i non-segmented inner door, handle right side, 100-127 V/50-60 Hz; Full stacking kit included	EP87-881
2 x CellXpert® C170i, inner door with 4 door segments, handle right side, 1 – 20 % oxygen control, water and humidity monitoring package, 100 – 127 V/50 – 60 Hz; Set of copper water tray and 4 copper shelves; Stacking stand, upper & lower frame	EP87-905
2 x CellXpert® C170i, non-segmented inner door, handle left side, 100-127 V/50-60 Hz; Full stacking kit included.	EP87-882



Make a qualifying CO₂ Incubator or Multipurpose Centrifuge and receive a FREE Cell Culture Consumables Box!

Free Consumables Box Includes:

- › 2x Eppendorf Serological Pipets, 5 mL, 4x100
- › 2x Eppendorf Serological Pipets, 10 mL, 4x100
- › 2x Eppendorf Serological Pipets, 25 mL, 4x50
- › Easypet® 3, 1-channel, 0.1-100 mL
- › 2x Eppendorf Conical Tubes Starter Pack, 25 mL
- › 2x Eppendorf Conical Tubes, SnapTec® cap, 25 mL

Contact to your Eppendorf Account Manager to take advantage today!

[eppendorf.com/us-en/salesreps](https://www.eppendorf.com/us-en/salesreps)



Eppendorf North America, Inc.
Phone: 800-645-3050

www.eppendorf.com